

[See all 9 Products in Family](#)

## 5X Nikon CF IC Epi Plan TI Michelson Interferometry Objective

See More by [Nikon](#)



Nikon Interferometry Objectives

Stock #59-311 [CONTACT US](#)

1 [-](#) [+](#) C\$7,133.<sup>00</sup>

[ADD TO CART](#)

### Volume Pricing

Qty 1+	C\$7,133.00 each
Need More?	<a href="#">Request Quote</a>

### Product Downloads

Field of View, 20 Diameter Field Eyepiece (mm):  
4.00

### General

Model Number:  
MUL42051

Compatible Tube Lens Focal Length (mm):  
Focal Length: 200mm

Type:  
Microscope Objective

Style:

Infinity Corrected

**Manufacturer:**  
Nikon

**Note:**  
Side port hosts the reference mirror and adjustment knob

## Physical & Mechanical Properties

**Field of View (mm):**  
4.00

**Length excluding Threads (mm):**  
39.70

**Maximum Diameter (mm):**  
31

**Weight (g):**  
280.00

## Optical Properties

**Horizontal Field of View, 1/2" Sensor:**  
1.28mm

**Horizontal Field of View, 2/3" Sensor:**  
1.76mm

**Focal Length FL (mm):**  
40.00

**Magnification:**  
5X

**Numerical Aperture NA:**  
0.13

**Resolving Power ( $\mu\text{m}$ ):**  
2.1

**Depth of Field ( $\mu\text{m}$ ):**  
16.20

**Field of View, 25 Diameter Field Eyepiece (mm):**  
5.00

**Working Distance (mm):**  
9.3

**Field Number (mm):**  
25

**Parfocal Length (mm):**  
49

**Immersion Liquid:**  
N/A

## Sensor

**Maximum Sensor Format:**  
2/3"

## Threading & Mounting

**Mounting Threads:**  
RMS / 20.32mm x 36 TPI

## Regulatory Compliance

**RoHS 2015:**  
[Compliant](#)

**Certificate of Conformance:**  
[View](#)

## Product Details

- Suitable for Non-Contact Optical Profiling
- Mchelson and Mrau Objectives Available
- Infinity Corrected 200mm Tube Lens Design

Nikon Interferometry Objectives are used in non-contact optical profile measurement devices to obtain surface maps and surface measurement parameters. They can be used to examine surface topography with very high precision—to within a fraction of the wavelength of light. In these objectives, a light beam passes through a beamsplitter, which directs the light to both the surface of the sample and a built-in reference mirror. The light reflected from these surfaces recombines and a fringe interference pattern is formed. Nikon Interferometry Objectives are available in magnifications ranging from 2.5X to 100X.

The Mchelson objectives provide comparatively longer working distances, wider fields of view and larger depth of field, whereas the Mrau objectives are used in applications requiring higher magnification and/or numerical apertures. Use of the Nikon 200mm tube lens ([#58-520](#)) allows these objectives to be integrated to a C-Mount camera.