

[See all 4 Products in Family](#)

# Olympus SZ61 Stereo Microscope Boom Stand System

See More by [Olympus](#)



SZ61 Stereo Microscope (eyepiece, objective, and stand each sold separately)

Stock **#89-352** [CONTACT US](#)

1  **C\$5,950<sup>00</sup>**

**ADD TO CART**

Volume Pricing	
Qty 1-4	<b>C\$5,950.00</b> each
Qty 5+	<b>C\$5,649.00</b> each
Need More?	<a href="#">Request Quote</a>

Product Downloads

## General

**Manufacturer:**

Olympus

## Regulatory Compliance

**Certificate of Conformance:**

[View](#)

## Product Details

- Simplify Component Selection for [SZ61](#) and [SZX7](#) Microscopes
- Systems Ready to use with Camera of your Choice
- Boom Stand Systems Include Ball Bearing Boom Stand

Preconfigured Olympus Stereo Microscope Systems include an Olympus SZ61 or SZX7 microscope and some of the most common components used to complete a microscope system. The [Olympus SZ61 Microscope](#) utilizes a Greenough optical system for improved image quality and superior field flatness. The [Olympus SZX7 Microscope](#) features video compatibility and a Galilean optical system for high resolution imaging. Boom stand systems include a high stability Ball Bearing Room Stand ([#54-121](#)). Microscope System includes a basic reflected light stand and is ready to be used with a [camera](#) of your choice. Add a camera (such as [Microscopy Camera](#)) to complete the video microscope system\*. See the Technical Information tab for parts that are included in each configuration.

\***Note:** Additional connection cables and/or external power supply may be required.

## Technical Information

Olympus SZ61 Stereo Microscope Boom Stand System ( <a href="#">#89-352</a> )	
Description	Stock No.
SZ61 Stereo Microscope Body with 45° Trinocular (C-Mount)	<a href="#">#88-126</a>
10X Widefield Eyepieces	<a href="#">#88-108</a>
0.3X Objective	<a href="#">#88-127</a>
Rack & Pinion Focusing Mount	<a href="#">#33-195</a>
32mm Fixed Column for Boom Stand	<a href="#">#55-618</a>

Olympus SZ61 Stereo Microscope System ( <a href="#">#89-354</a> )	
Description	Stock No.
SZ61 Stereo Microscope Body with 45° Trinocular (C-Mount)	<a href="#">#88-126</a>
10X Widefield Eyepieces	<a href="#">#88-108</a>
0.3X Objective	<a href="#">#88-127</a>
Basic Reflected Light Stand	<a href="#">#88-207</a>
Your choice of <a href="#">Camera</a> and <a href="#">Lighting</a> *	Not included*

Olympus SZX7 Stereo Microscope Boom Stand System ( <a href="#">#89-353</a> )	
Description	Stock No.
SZX7 Stereo Microscope Body	<a href="#">#88-103</a>
30° Binocular Observation Tube	<a href="#">#88-106</a>
10X Widefield Eyepieces	<a href="#">#88-108</a>
1X Achromatic Objective Lens	<a href="#">#88-112</a>
Rack & Pinion Focusing Mount	<a href="#">#33-195</a>
32mm Fixed Column for Boom Stand	<a href="#">#55-618</a>

Olympus SZX7 Stereo Microscope System ( <a href="#">#89-355</a> )	
Description	Stock No.
SZX7 Stereo Microscope Body	<a href="#">#88-103</a>
30° Trinocular Observation Tube	<a href="#">#88-107</a>
10X Widefield Eyepieces	<a href="#">#88-108</a>
1X Achromatic Objective Lens	<a href="#">#88-112</a>
1X Camera adapter with C-Mount Coupler	<a href="#">#88-209</a>
Basic Reflected Light Stand	<a href="#">#88-207</a>
Your choice of <a href="#">Camera</a> and <a href="#">Lighting</a> *	Not included*

\*Choose your preferred Camera (such as [Microscopy Camera](#)), additional Lighting (such as [Microscopy Lighting](#)) and other optional components from our wide selection to complete a video microscope system for industrial inspection or life science application. These are sold separately. Additional connection cables and/or external power supply may be required

