

[See all 4 Products in Family](#)

# Olympus SZX7 Stereo Microscope System

See More by [Olympus](#)



SZX7 Zoom Stereo Microscope - accessories shown not included

Stock **#89-355** [CONTACT US](#)

-

1

+

C\$10,430<sup>.00</sup>

ADD TO CART

| Volume Pricing |                               |
|----------------|-------------------------------|
| Qty 1-4        | C\$10,430.00 each             |
| Qty 5+         | C\$9,905.00 each              |
| Need More?     | <a href="#">Request Quote</a> |

## Product Downloads

### General

Manufacturer:

Olympus

### Regulatory Compliance

Certificate of Conformance:

[View](#)

Product Details

- Simplify Component Selection for [SZ61](#) and [SZX7](#) Microscopes
- Systems Ready to use with Camera of your Choice
- Boom Stand Systems Include Ball Bearing Boom Stand

Preconfigured Olympus Stereo Microscope Systems include an Olympus SZ61 or SZX7 microscope and some of the most common components used to complete a microscope system. The [Olympus SZ61 Microscope](#) utilizes a Greenough optical system for improved image quality and superior field flatness. The [Olympus SZX7 Microscope](#) features video compatibility and a Galilean optical system for high resolution imaging. Boom stand systems include a high stability Ball Bearing Room Stand ([#54-121](#)). Microscope System includes a basic reflected light stand and is ready to be used with a [camera](#) of your choice. Add a camera (such as [Microscopy Camera](#)) to complete the video microscope system\*. See the Technical Information tab for parts that are included in each configuration.

\***Note:** Additional connection cables and/or external power supply may be required.

Technical Information

|                                                                              |                         |
|------------------------------------------------------------------------------|-------------------------|
| Olympus SZ61 Stereo Microscope Boom Stand System ( <a href="#">#89-352</a> ) |                         |
| Description                                                                  | Stock No.               |
| SZ61 Stereo Microscope Body with 45° Trinocular (C-Mount)                    | <a href="#">#88-126</a> |
| 10XWidefield Eyepieces                                                       | <a href="#">#88-108</a> |
| 0.3X Objective                                                               | <a href="#">#88-127</a> |
| Rack & Pinion Focusing Mount                                                 | <a href="#">#33-195</a> |
| 32mm Fixed Column for Boom Stand                                             | <a href="#">#55-618</a> |

|                                                                      |                         |
|----------------------------------------------------------------------|-------------------------|
| Olympus SZ61 Stereo Microscope System ( <a href="#">#89-354</a> )    |                         |
| Description                                                          | Stock No.               |
| SZ61 Stereo Microscope Body with 45° Trinocular (C-Mount)            | <a href="#">#88-126</a> |
| 10XWidefield Eyepieces                                               | <a href="#">#88-108</a> |
| 0.3X Objective                                                       | <a href="#">#88-127</a> |
| Basic Reflected Light Stand                                          | <a href="#">#88-207</a> |
| Your choice of <a href="#">Camera</a> and <a href="#">Lighting</a> * | Not included*           |

|                                                                              |                         |
|------------------------------------------------------------------------------|-------------------------|
| Olympus SZX7 Stereo Microscope Boom Stand System ( <a href="#">#89-353</a> ) |                         |
| Description                                                                  | Stock No.               |
| SZX7 Stereo Microscope Body                                                  | <a href="#">#88-103</a> |
| 30° Binocular Observation Tube                                               | <a href="#">#88-106</a> |
| 10XWidefield Eyepieces                                                       | <a href="#">#88-108</a> |
| 1XAchromatic Objective Lens                                                  | <a href="#">#88-112</a> |
| Rack & Pinion Focusing Mount                                                 | <a href="#">#33-195</a> |
| 32mm Fixed Column for Boom Stand                                             | <a href="#">#55-618</a> |

|                                                                      |                         |
|----------------------------------------------------------------------|-------------------------|
| Olympus SZX7 Stereo Microscope System ( <a href="#">#89-355</a> )    |                         |
| Description                                                          | Stock No.               |
| SZX7 Stereo Microscope Body                                          | <a href="#">#88-103</a> |
| 30° Trinocular Observation Tube                                      | <a href="#">#88-107</a> |
| 10XWidefield Eyepieces                                               | <a href="#">#88-108</a> |
| 1XAchromatic Objective Lens                                          | <a href="#">#88-112</a> |
| 1XCamera adapter with C-Mount Coupler                                | <a href="#">#88-209</a> |
| Basic Reflected Light Stand                                          | <a href="#">#88-207</a> |
| Your choice of <a href="#">Camera</a> and <a href="#">Lighting</a> * | Not included*           |

\*Choose your preferred Camera (such as [Microscopy Camera](#)), additional Lighting (such as [Microscopy Lighting](#)) and other optional components from our wide selection to complete a video microscope system for industrial inspection or life science application. These are sold separately. Additional connection cables and/or external power supply may be required

