

[See all 13 Products in Family](#)

## 10X Objective, CFI Plan Achromatic

See More by [Nikon](#)



Stock #75-353 NEW CONTACT US

-
1
+
C\$439<sup>00</sup>

ADD TO CART

Volume Pricing	
Qty 1+	C\$439.60 each
Need More?	<a href="#">Request Quote</a>

Product Downloads

**General**

**Model Number:**  
MRL00102

**Compatible Tube Lens Focal Length (mm):**  
Focal Length: 200mm

**Type:**  
Microscope Objective

**Style:**  
Infinity Corrected

**Manufacturer:**  
Nikon

## Physical & Mechanical Properties

2.2 **Field of View (mm):**

49.40 **Length excluding Threads (mm):**

27 **Maximum Diameter (mm):**

90 **Weight (g):**

## Optical Properties

N/A **Compatible Cover Glass Thickness (mm):**

0.64 **Horizontal Field of View, 1/2" Sensor:**

0.88 **Horizontal Field of View, 2/3" Sensor:**

10X **Magnification:**

0.25 **Numerical Aperture NA:**

10.5 **Working Distance (mm):**

22 **Field Number (mm):**

59.9 **Parfocal Length (mm):**

N/A **Immersion Liquid:**

## Sensor

2/3" **Maximum Sensor Format:**

## Threading & Mounting

M25 x 0.75 **Mounting Threads:**

## Regulatory Compliance

[View](#) **Certificate of Conformance:**

## Product Details

- Exceptional Flat-Field Imaging
- High Numerical Apertures and Oil Immersion Options Available
- Wide Magnification Range (1X to 100X)

Nikon's CFI Plan Achromat Objectives deliver exceptional flat-field imaging ensuring sharp, distortion-free clarity across the entire field of view, making these objectives ideal for both visual inspection and high-precision digital imaging. With high numerical apertures and specialized oil immersion options available, these objectives offer enhanced resolution and light-gathering capability for demanding high-magnification applications. Nikon CFI Plan Achromat Objectives are available in 1X up to 100X magnification, providing solutions for low-magnification overviews or detailed high-resolution imaging. Color corrected for the entire visible spectrum; these objectives are suitable for brightfield and fluorescence observation in routine lab work and photomicrography.

## Technical Information

