

TECHSPEC[®] Cy5.5™ Filter Cube Set - Nikon



Olympus Filter Cube #86-833



Stock **#67-011-NKN** 1 In Stock

-

1

+

C\$1,773^{.38}

ADD TO CART

Volume Pricing	
Qty 1-9	C\$1,773.38 each
Qty 10-25	C\$1,613.11 each
Qty 26-49	C\$1,522.58 each
Need More?	Request Quote

Product Downloads

SPECIFICATIONS

General

Type:

Fluorescence Filter Kit	
CY5.5	Compatible Fluorophore:
#67-037	Excitation Filter:
#67-039	Emission Filter:
#67-085	Dichroic Filter:
Nikon	Compatible Microscope Brand:
Optical Properties	
Fused Silica (Corning 7980)	Substrate: ☐
Hard Coated	Coating:
635 - 675	Excitation Wavelength (nm):
696 - 736	Emission Wavelength (nm):
685.00	Dichroic Cut-On Wavelength (nm):
635 - 736	Wavelength Range (nm):
Regulatory Compliance	
View	Certificate of Conformance:

PRODUCT DETAILS

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).