

TECHSPEC® Mounting Clamp for 2 - 8X Research-Grade Variable Beam Expander



Mbunting Clamp for 2 - 8X Research-Grade Variable Beam Expander, #89-295

Stock **#89-295** ☐ In Stock

- 1 + **C\$229.60**

ADD TO CART

Qty 1-4

C\$229.60

Qty 5+

C\$203.00

Volume Pricing

[Request Quote](#)

Get Product Downloads



Type:
Mounting Clamp

Regulatory Compliance

RoHS:
[Compliant](#)

