

TECHSPEC[®] SpectrumGreen Filter Cube Set - Nikon



Olympus Filter Cube #86-833



Stock **#86-373-NKN** [CONTACT US](#)

-

1

+

C\$1,890^{.00}

ADD TO CART

Volume Pricing	
Qty 1-9	C\$1,890.00 each
Qty 10+	C\$1,701.00 each
Need More?	Request Quote

Product Downloads

SPECIFICATIONS

General

Fluorescence Filter Kit

Type:

SpectrumGreen	Compatible Fluorophore:
#84-096	Excitation Filter:
#84-098	Emission Filter:
#67-080	Dichroic Filter:
Nikon	Compatible Microscope Brand:
Optical Properties	
Fused Silica (Corning 7980)	Substrate:
Hard Coated	Coating:
484 - 504	Excitation Wavelength (nm):
517 - 537	Emission Wavelength (nm):
506.00	Dichroic Cut-On Wavelength (nm):
484 - 537	Wavelength Range (nm):
Regulatory Compliance	
View	Certificate of Conformance:

PRODUCT DETAILS

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).