

[See all 15 Products in Family](#)

## Standard In-Line Attachment (Optimized at 2X-10X), InfiniTube

See More by [Infinity Photo-Optical Company](#)



Standard In-Line Attachment (Optimized at 2X-10X), InfiniTube

Stock **#56-126** **1 In Stock**

⊖ 1 ⊕ C\$1,736<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1+	C\$1,736.00 each
Need More?	<a href="#">Request Quote</a>

### Product Downloads

### General

Collector I **Model Number:**

Fixed Magnification Lens **Type:**

### Optical Properties

Objective Required **Primary Magnification PMAG:**

## Regulatory Compliance

Compliant

RoHS 2015:

Compliant

Reach 224:

[View](#)


[Certificate of Conformance:](#)

## Product Details


- Enable Dynamic Optical Focusing and Reduce Spherical Aberration
- Compatible With All Infinity Corrected Objectives
- Color Matched for Major Brand Objectives

InFocus™ Dynamic Optical Focusing Systems are used with microscopes to enable dynamic optical focusing, and drastically reduce the spherical aberration caused by the sample's cover glass thickness or suspension fluid. These dynamic optical focusing systems eliminate the need for correction collars and expensive motorized computer activation. InFocus™ Dynamic Optical Focusing Systems allows micromanipulation of the microscope, and can be integrated with any infinity corrected objectives. Variations in magnification are minimal over the range of focal translation. InFocus™ Dynamic Optical Focusing System modules are color corrected for Mitutoyo, Achrovid, Zeiss, Leica, Nikon, or Olympus objectives and can be connected to most major brand microscopes using a series of stand adapters (sold separately).

## Technical Information



### InfiniTube™ In-Line Attachment



Units: mm