

[See all 5 Products in Family](#)

**TECHSPEC® T-Mount to Olympus (RMS/DIN) Adapter**



T-Mount Objective Ring Mount, #52-301

Stock **#52-301** **14 In Stock**

⊖ 1 ⊕ C\$83<sup>30</sup>

**ADD TO CART**

Volume Pricing	
Qty 1-4	C\$83.30 each
Qty 5-9	C\$74.06 each
Qty 10-25	C\$68.74 each
Need More?	<a href="#">Request Quote</a>

Product Downloads

Objective Mount	Type:
48.0	Outer Diameter (mm):
RMS/DIN (0.8" - 36 TPI)	Objective Thread:
T-Mount	Thread Type:
	Length (mm):

11.5      **Extension Length (mm):**

Objective Adapter      **Function:**

Objective      **Type of Optics:**

**Regulatory Compliance**

**Compliant**      **RoHS 2015:**

**View**      **Certificate of Conformance:**

**Compliant**      **REACH 241:**


**Product Details**

- Integrate Microscope Objectives to Create Digital Microscope Systems
- Feature Male C-Mount, T-Mount, or ID1 (SM1) Thread and Female Objective Thread
- Black Anodized Aluminum Construction

TECHSPEC® Tube System Objective Adapters integrate microscope objectives into C-Mount, T-Mount, or ID1 (SM1) systems. These adapters offer female threads for multiple standard objective thread types, including Olympus (RMS/DIN), Mitutoyo (M26 x 0.706), and Nikon (M25 x 0.75). A black anodized finish increases these adapters' durability while reducing stray light and reflections within the system. TECHSPEC® Tube System Objective Adapters allow for the personalization and customization of imaging systems. They can be used to construct digital microscope systems by allowing microscope objectives to be integrated into digital camera systems.

**Technical Information**

**Configuration Example - C-Mount Digital Microscope System with Infinity Corrected Objectives**



Dimensions			
Compatible with	Overall Size (Dia. x L)	Required Tube Length (A)	
Mitutoyo MT-1/MT-2	50 x 42mm	190mm	
Mitutoyo MT-L/MT-L4	50 x 40.4mm	150mm	
Mitutoyo MT-4/MT-40	30 x 152.5mm	0mm	
Olympus WF (#36-401)	56 x 47.8mm	123mm	
Nikon (#58-520)	50 x 37.7mm	130mm	